

CT Bony Pelvis without contrast

Clinical Indication	Comment
Trauma, pre-operative planning, pain	

ACQUISITION

Contrast Phase	
Extent	
Scan direction	
Acquisition technique	
IV Contrast technique	NA
IV Contrast rate	NA
Saline flush	
Delay	
SmartPrep/ Bolus Pro	
Patient Age	All ages
Patient Size	All sizes
Dose Index Range	Action level for investigation: >500 mGy total CTDIvol, >5000 mGy*cm total DLP

RECONSTRUCTION

Source	Name of series	Reconstruction/Reformat
Raw data	SAG	Sagittal
Raw data	COR	Coronal
Raw data	3D	3D Volume rendered

Keywords

Trauma
pre-operative planning
pain